



ENDOSCOPY SERVICES

Competitive prices in today's economy to meet your needs

**Flexible Video Laparoscopes EndoEYE Deflectable Tip
LTF-VH High Definition and LTF-VP, 5mm, 0 Degree:**

Deflecting Tip Repair Level 1

- Angulations Adjusted
- Control Lever Replaced
- Bending Rubber Replaced
- Cover Lens Polished
- Electrical Connector Pins Cleaned
- Light Guide Prong Cleaned
- Moisture Removed
- Complete Switch Assembly Repaired, including Ribbon and Buttons

Deflecting Tip Repair Level 2

- Electrical Circuit Board Repaired
- Light Guide/ Electrical Connector Tube Replaced
- Light Guide Connector Prong End Lens Cap Repaired/Replaced
- Angulations Adjusted, including Stoppers

Deflecting Tip Repair Level 3

- Control Body Prong End Replacement
- Light Guide Bundle cleaned and refitted/replaced
- Electrical Circuit Board Replaced-with entire components
- Angulation System Overhaul
- Electrical Connector Overhaul
- Light Guide Prong Adapter repaired/ replaced

Deflecting Tip Repair Level 4

- CCD Chip re-wire
- Optical System Repaired/Replaced
- CCD Chip Replacement- all components, including optical system assembly

*Limitations only possible in the rare event that a CCD can't be repaired and replacement can't be found in a timely manner.

**90 days warranty 1-2 days minor repairs / 3-5 days major repairs

Thank you for choosing Elite Endoscopy Services for all your repairs.

Office: 954-543-4333

Fax: 954-543-4331

info@eeservicescorp.com

www.ELITEendoscopyservices.com