



ENDOSCOPY SERVICES

Competitive prices in today's economy to meet your needs

## ENDOYEYE RIGID REPAIR LEVELS

Elite Endoscopy Services will meet all your endoscope needs for all medical fields. Our technicians are highly trained and skilled to repair your endoscopes.

### **Rigid Video Laparoscopes and HD Endoeye Telescopes:**

#### **Repair Level 1: Rigid Video Laparoscopes Endoeye**

- o Electrical connector repaired
- o Switch covers repaired
- o Bent shaft repaired
- o Minor dents and scratches on shaft removed
- o Lens cover polished
- o Light guide lens prong replacement
- o Cuts and/ or small holes in light guide tube repaired

#### **Repair Level 2: Rigid Video Laparoscopes Endoeye**

- o Electrical circuit board repair
- o Light guide tube/ Electrical connector tube replacement
- o Switch complete assembly repair/ replacement
- o Wire harness re-wired/ replacement
- o Electrical housing paddle replacement
- o Inner and Outer tube replacement (re-fiber) with new desiccant installed
- o Insertion tube shaft window at distal tip end replacement

#### **Repair Level 3: Rigid Video Laparoscopes Endoeye**

- o CCD Chip re-wire with new components installed
- o Electrical Circuit Board replacement
- o Optical System repair/ replacement
- o CCD chip negative lens replacement

#### **Repair Level 4: Rigid Video Laparoscopes Endoeye**

- o CCD Chip complete overhaul with optical system replacement

\*Limitations only possible in the rare event that a CCD can't be repaired and replacement can't be found in a timely manner.

\*\*90 days warranty 1-2 days minor repairs / 3-5 days major repairs

**Thank you for choosing Elite Endoscopy Services for all your repairs.**

Office: 954-543-4333

Fax: 954-543-4331

[info@eeservicescorp.com](mailto:info@eeservicescorp.com)

[www.ELITEendoscopyservices.com](http://www.ELITEendoscopyservices.com)